

ABCB1 GENE VARIANT C1236T (MDR1)

ORDERING INFORMATIONS

REF: FGC-003-25
 RDM Code: 1875566/R
 CND Code: W0106010499
 Tests: 25
 Reactions: 31
 Manufacturer: BioMol Laboratories s.r.l.

CONTENTS OF THE KIT

The kit consists of: reagents for Real-Time PCR amplification
 *the reagents for the extraction of genomic DNA are not supplied in the kit.

For in vitro diagnostic use



PRODUCT CHARACTERISTICS

The FGC-003 kit allows the characterization of the genetic variant C1236T of the ABCB1 gene (rs1128503) by amplification with oligonucleotides and specific probes (allele-specific genotyping) and subsequent detection with qPCR-Real-time. Kit optimized for Real-Time PCR instrumentation Biorad CFX96, Biorad Opus Dx, Agilent AriaDx.

SCIENTIFIC BACKGROUND

Docetaxel and paclitaxel are cytotoxic taxanes that inhibit mitosis causing the death of cancer cells. They are mainly used in the treatment of breast, ovarian and lung cancer. For taxanes, the ABCB1 gene is considered one of the best candidates to explain variations in clinical responses and toxicity.

Doxorubicin, an anthracycline widely used as mono- or combination therapy in the treatment of solid tumors including breast cancer, is also the substrate of P-gp. Significantly impaired clearance and decreased plasma concentration of doxorubicin have been observed in patients with one of the three above described polymorphisms of the ABCB1 gene. Irinotecan, a topoisomerase I inhibitor, plays a significant role in the treatment of colorectal cancer in monotherapy or in combination with 5-FU. The antitumor activity is mainly due to the SN-38 metabolite which is metabolized to SN-38G, which has 1/100 of the antitumor activity and is practically inactive.

CLINICAL SIGNIFICANCE

Pharmacogenetic screening and/or drug-specific phenotyping of cancer patients eligible for treatment with chemotherapy drugs can identify patients likely to be reactive or resistant to proposed drugs. Similarly, identification of patients with an increased risk of developing toxicity allows for dose adaptation or application of other targeted therapies. Polymorphisms in genes encoding drug efflux transporters, such as P-glycoprotein, can affect the absorption and excretion of anticancer drugs. It is encoded by the multidrug resistance 1 (MDR1) gene (ABCB1, ATP-binding cassette transporter superfamily B member 1) located on chromosome 7q21. There are three main polymorphisms affecting P-gp activity: the c.2677G>T/A polymorphism in exon 21 which causes a substitution in the amino acid sequence Ala (G)/Ser (T) or Thr (A), with consequent possible increase in enzyme function. The second polymorphism is in exon 26 at position c.3435C>T, resulting in more than twofold expression of P-gp. The third C1236T polymorphism in exon 12 does not directly affect P-gp expression but has an indirect effect as it alters the stability of the mRNA encoding the protein.

§ Clinical utility of ABCB1 genotyping for preventing toxicity in treatment with irinotecan. *Pharmacol Res.* 2018 Oct; 136:133-139. doi:10.1016/j.phrs.2018.08.026.

Epub 2018 Sep 11.

§ Genotypes Affecting the Pharmacokinetics of Anticancer Drugs. *Clin Pharmacokinet.* 2017, Apr; 56(4):317-337. doi: 10.1007/s40262-016-0450-z. Review.

§ Influence of the ABCB1 polymorphisms on the response to Taxane-containing chemotherapy: a systematic review and meta-analysis. *Cancer Chemother Pharmacol.* 2018, Feb; 81(2):315-323. doi: 10.1007/s00280-017-3496-1. Epub 2017 Dec 5.

ABCB1 GENE VARIANT C1236T (MDR1)

ORDERING INFORMATIONS

REF: FGC-003-25
 RDM Code: 1875566/R
 CND Code: W0106010499
 Tests: 25
 Reactions: 31
 Manufacturer: BioMol Laboratories s.r.l.

CONTENTS OF THE KIT

The kit consists of: reagents for Real-Time PCR amplification
 *the reagents for the extraction of genomic DNA are not supplied in the kit.

For in vitro diagnostic use



CONTENTS OF THE KIT

DESCRIPTION	LABEL	VOLUME	STORAGE
		FGC-003-25	
Mix oligonucleotides and probes	Mix 10X C1236T ABCB1	1 x 85 µl	-20°C
Mix buffer and Taq-polymerase enzyme	Mix Real-Time PCR 2X	1 x 425 µl	-20°C
Deionized H ₂ O	Deionized H ₂ O	2 x 1 ml	-20°C
Genomic DNA or recombinant DNA	Control +1	1 x 22 µl	-20°C
Genomic DNA or recombinant DNA	Control +2	1 x 22 µl	-20°C
Genomic DNA or recombinant DNA	Control +2	1 x 22 µl	-20°C

TECHNICAL CHARACTERISTICS

COD. FGC-003-25

STABILITY	18 months
REAGENTS STATUS	Ready to use
BIOLOGICAL MATRIX	Genomic DNA extracted from whole blood, tissue, cells
POSITIVE CONTROL	Recombinant DNA for at least 3 analytical sessions
TECHNOLOGY	Real-time PCR; oligonucleotides and specific probes; 2 FAM/HEX fluorescence channels
VALIDATED INSTRUMENTS	Biorad CFX96 Dx, Biorad Opus Dx e Agilent AriaDx
RUNNING TIME	85 min
THERMAL CYCLING PROFILE	1 cycle at 95 °C (10 min); 45 cycles at 95 °C (15 sec) + 60 °C (60 sec)
ANALYTICAL SPECIFICITY	Absence of non-specific pairings of oligonucleotides and probes; absence of cross-reactivity
ANALYTICAL SENSITIVITY : LIMIT OF DETECTION (LOD)	≥ 0,016 ng of DNA
ANALYTICAL SENSITIVITY : LIMIT OF BLANK (LOB)	0% NCN
REPRODUCIBILITY	99,9%
DIAGNOSTIC SPECIFICITY / DIAGNOSTIC SENSITIVITY	100%/98%