





# SCHEMA TECNICA

## 8-Channel Adjustable Volume Pipette



Multi-Channel Adjustable Volume Pipette

**Precipette™**



# SCHEMA TECNICA

## 8-Channel Adjustable Volume Pipette

### Color-coded control button



0.5ul-300ul  
Wide volume range

Autoclavable lower section  
at **121°C**



### Ergonomic design

Comfortable grip  
Fatigue-free pipetting



### Easy to clean and maintain



| Model No | Volume Range | Increment | Testing volume |
|----------|--------------|-----------|----------------|
| MP10     | 0.5-10 µl    | 0.01 µl   | 10, 5, 1       |
| MP50     | 5-50 µl      | 0.05 µl   | 50, 25, 5      |
| MP300    | 30-300 µl    | 0.5 µl    | 300, 150, 30   |