





SCHEMA TECNICA


Fully Autoclavable Adjustable Volume Pipette





SCHEMA TECNICA

Fully Autoclavable Adjustable Volume Pipette



Model No	Volume Range	Increment	Test volume
PPT0101	0.1-2.5 μ l	0.002 μ l	0.25, 1.25, 2.5 μ l
PPT0102	0.5-10 μ l	0.01 μ l	1, 5, 10 μ l
PPT0103	2-20 μ l	0.02 μ l	2, 10, 20 μ l
PPT0104	5-50 μ l	0.05 μ l	5, 25, 50 μ l
PPT0105	10-100 μ l	0.1 μ l	10, 50, 100 μ l
PPT0106	20-200 μ l	0.2 μ l	20, 100, 200 μ l
PPT0107	100-1000 μ l	1 μ l	100, 500, 1000 μ l
PPT0108	0.5-5 ml	5 μ l	500, 2500, 5000 μ l
PPT0109	1-10 ml	10 μ l	1, 5, 10 ml